


*Hi, I am Aurocoat Plus.*



**OPHTHALMIC VISCOSURGICAL DEVICE**  
(Sodium hyaluronate 3%w/v with Chondroitin Sulphate 4%w/v)



## What am I?

Aurocoat plus ophthalmic viscosurgical device is a sterile, isotonic, pyrogen-free, high viscosity, medium molecular weight of highly purified grade of sodium hyaluronate, and sodium chondroitin sulphate dissolved in a buffered physiological solution for intraocular application. It is a clear solution in a ready-to-use syringe for single use. It is packed in a sterile blister or heat sealed pouch.

## Properties

Composition		
Sodium hyaluronate	:	30mg/ml (3%w/v)
Chondroitin sulphate	:	40mg/ml (4%w/v)
Molecular weight		
Sodium hyaluronate	:	1M daltons
Chondroitin sulphate	:	> 22000 daltons
pH	:	6.5 - 7.8
Osmolality	:	250 - 400 mOsm/Kg
Viscosity @ 250 C	:	20000 - 60000 cPs
Refractive index @ 250 C	:	1.33 - 1.35
Specific gravity	:	1.000 - 1.040
Spectral transmittance	:	>90%

## What I do?

*During cataract surgery*

- I protect the corneal endothelium
- I provide anterior chamber stability
- I aid in superior post operative results
- I enter with ease and leave as you want

## Supply

*Sterile 1ml pre-filled glass syringes*

## Storage

*Store between 2° and 8°C protected from light. Do not freeze.*

Information published in this catalogue is subject to change without notification



No.1, Sivagangai Main Road, Veerapanjan, Madurai 625 020, India  
Phone: +91 452 2446200 Fax: +91 452 2446200  
Email: info@aurolab.com Web:www.aurolab.com

